




***Kehr Surgical Private Limited***



**Kehr Surgical Private Limited**  
*The smart choice for sharp blades*



**KiaTo**

### Our company

Kehr Surgical Private Limited has been manufacturing superior quality surgical blades since 1970. Exports have contributed significantly to sales since 1975. We export to more than 80 countries.

The confidence of leading Surgeons and hospitals worldwide is testimony to the Company's ability to supply high quality blades at very competitive prices.

We are proud of our ability and history of not only producing high quality world standard surgical blades but also having the distinction of producing a wide range of specialized blades.

### Our products



CE Marked surgical blades and disposable scalpels are available in both carbon steel and stainless steel. The blades are manufactured in the following sizes: 10,10A, 11,12, 15, 18, 20, 21, 22, 23, 24, 25 and 36. The Company also has the ability to manufacture other sizes as per customer requirements. The blades are packaged in individual aluminum foil pouch, which are then packed in lots of one hundred and sterilized by gamma radiation. The scalpels pouches are packed in boxes of ten.

We also offer Stich Cutters, Safety Scalpel, Post Mortem Blades, Gouge Blades, Podiatry blades and special purpose blades

### Our manufacturing

At Kehr Surgical Private Limited product quality and process parameters are stringently monitored and controlled at each and every stage of production. Bar codes is used for identification and product traceability.

The facility is ISO 13485:2003/ ISO 13485:2012 certified by DNV GL Nemko. We recognize the importance of continued improvements in process capability in order to achieve the objective of keeping up our reputation as a leading quality supplier. The medical devices are packed in a class 1000 clean room.



C34 Panki Industrial Estate  
Kanpur 208 022 UP India  
Tel: 91-512-2692507 2692508 Fax: 91-512-2692509  
Email: info@kehrsurgical.com  
Website: www.kehrsurgical.com



EDI TI QN 04/ Dec 17 KJ ATO