




Kehr Surgical Private Limited



Kehr Surgical Private Limited
The smart choice for sharp blades



KIATO

Our company

Kehr Surgical Private Limited has been manufacturing superior quality surgical blades since 1970. Exports have contributed significantly to sales since 1975. We export to more than 80 countries.

The confidence of leading Surgeons and hospitals worldwide is testimony to the Company's ability to supply high quality blades at very competitive prices.

We are proud of our ability and history of not only producing high quality world standard surgical blades but also having the distinction of producing a wide range of specialized blades.

Our products



CE Marked surgical blades and disposable scalpels are available in both carbon steel and stainless steel. The blades are manufactured in the following sizes: 10,10A, 11,12, 15, 18, 20, 21, 22, 23, 24, 25 and 36. The Company also has the ability to manufacture other sizes as per customer requirements. The blades are packaged in individual aluminum foil pouch, which are then packed in lots of one hundred and sterilized by gamma radiation. The scalpels pouches are packed in boxes of ten.

We also offer Stich Cutters, Safety Scalpel, Post Mortem Blades, Gouge Blades, Podiatry blades and special purpose blades

Our manufacturing

At Kehr Surgical Private Limited product quality and process parameters are stringently monitored and controlled at each and every stage of production. Bar codes is used for identification and product traceability.

The facility is ISO 13485:2016 certified by DNV GL Presafe AS. We recognize the importance of continued improvements in process capability in order to achieve the objective of keeping up our reputation as a leading quality supplier. The medical devices are packed in a class 1000 clean room.



C34 Panki Industrial Estate
Kanpur 208 022 UP India
Tel: 91-512-2692507 2692508 Fax: 91-512-2692509
Email: info@kehrsurgical.com
Website: www.kehrsurgical.com

Kiato Brand Carbon Steel Surgical Blade

Size	Reorder No:	Size	Reorder No:
09	11189	20	11120
10	11110	21	11121
10A	11109	22	11122
11	11111	22A	11130
11P	11161	23	11123
12	11112	24	11124
12B	11163	24D	11151
12D	11125	25	11128
13	11113	25A	11186
14	11114	26	11133
15	11115	27	11139
15A	11164	34	11134
15B	11165	36	11129
15C	11126	36D	11136
15T	11145	60	11166
16	11116	PM40	11140
17	11117	PM40B	11141
18	11127	PM60	11160
19	11119	PM60B	11162