



Skin Graft Blade – 158 mm

Compatible with Watson, Cobbett & Braithwaite Skin Graft Handles

Product Overview

The KIATO Skin Graft Blade (158 mm) is a sterile, single-use medical device designed for the harvesting of split-thickness skin grafts when used with compatible skin graft handles. Manufactured from medical-grade stainless steel and supplied sterile, the blade is produced under ISO 13485 certified quality management systems to support consistent and reliable clinical use.

Compliance & Use

This device is intended for use by trained healthcare professionals only. The blade is supplied sterile for single use and should not be reused or reprocessed. Certification details and documentation are available upon request for distributor qualification and tenders.

Product Specifications

Product Name	Skin Graft Blade
Brand	KIATO
Product Reference	12400
Intended Use	Harvesting of split-thickness skin grafts
Usage Type	Single Use / Disposable
Blade Material	Medical-grade Stainless Steel
Blade Length	158 mm
Blade Thickness	Nominal 0.25 mm
Blade Edge Design	Precision-ground, symmetric double-facet cutting edge
Handle Compatibility	Watson, Cobbett and Braithwaite Skin Graft Handles
Sterility	Supplied Sterile (Gamma Sterilised)
Manufacturing Standard	Manufactured under ISO 13485 certified quality management system
Packaging	Individually sterile packed in easy-peel packs
Shelf Life	5 years from date of manufacture
Export Markets	Global

Clinical Applications

Used by trained healthcare professionals for split-thickness skin graft harvesting in burn wound management, reconstructive and plastic surgery, trauma-related skin loss, chronic wound management, and post-oncological skin reconstruction. **Procedure selection and technique are determined by the clinician.**



Product Image & Dimensions



Illustrative representation for reference only.

Dimension	Specification
Blade Length	158 mm
Blade Thickness	Nominal 0.25 mm
Handle Compatibility	Watson, Cobbett & Braithwaite



Marketed By:



Kehr Surgical Private Limited
The smart choice for sharp blades

C-34 Panki Industrial Estate, Site 1,
Kanpur - 208 022, U.P. India
TEL : 0512-4058442 / 73111 82746
www.kehrsurgical.com
CIN No. U33111UP1970PTC003356
UP4842MD00023.



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland



aditya
DISPOMED

A subsidiary of Kehr Surgical Pvt. Ltd.



Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana I22050 India

Tel. +91 I24 47 64 900
Fax +91 I24 47 64 911

info@adityadispomed.com
www.adityadispomed.com