



Skin Graft Blade – 158 mm

Compatible with Watson, Cobbett & Braithwaite Skin Graft Handles

Product Overview

The SURGEON Skin Graft Blade (158 mm) is a sterile, single-use medical device designed for the harvesting of split-thickness skin grafts when used with compatible skin graft handles. Manufactured from medical-grade stainless steel and supplied sterile, the blade is produced under ISO 13485 certified quality management systems to support consistent and reliable clinical use.

Product Specifications

Product Name	Skin Graft Blade
Brand	Surgeon
Product Reference	12450
Intended Use	Harvesting of split-thickness skin grafts
Usage Type	Single Use / Disposable
Blade Material	Medical-grade Stainless Steel
Blade Length	158 mm
Blade Thickness	Nominal 0.25 mm
Blade Edge Design	Precision-ground, symmetric double-facet cutting edge
Handle Compatibility	Watson, Cobbett and Braithwaite Skin Graft Handles
Sterility	Supplied Sterile (Gamma Sterilised)
Manufacturing Standard	Manufactured under ISO 13485 certified quality management system
Packaging	10 Units per Box
Shelf Life	5 years from date of manufacture

Clinical Applications

Used by trained healthcare professionals for split-thickness skin graft harvesting in burn wound management, reconstructive and plastic surgery, trauma-related skin loss, chronic wound management, and post-oncological skin reconstruction. Procedure selection and technique are determined by the clinician.



Product Image & Dimensions



Illustrative representation for reference only.

Dimension	Specification
Blade Length	158 mm
Blade Thickness	Nominal 0.25 mm
Handle Compatibility	Watson, Cobbett & Braithwaite



Marketed By:



Kehr Surgical Private Limited
The smart choice for sharp blades
 C-34 Panki Industrial Estate. Site 1.
 Kanpur - 208 022, U.P. India
 Customer Care No. 0512-4058442 / 73111 82746
 Email: info@kehrsurgical.com
 CIN No. U33111UP1970PTC003356